



**GAWANDE COLLEGE OF
PHARMACY. Sakhar Kherda**
Sakharkherda, Lavala Road, Tq.
Sindkhed Raja
Dist. Buldhana



Department: Pharmaceutics

Name of the instrument: **Ball Mill.**

S.O.P. No.: GCP/PH/062


Stock Register No.: PH-00106/01

Effective date: 20/02/15

Review date: 27/02/15

Procedure:

- Ensure the mill is properly clean and dry.
- The drug is loaded in ball mill.
- Switch on the button.
- Adjust the speed of rotation.
- Maintain the speed of for specific time.
- Switch of the button.
- After specific time remove the drug and pass through sieve to get required drug size.



Ball Mill.

<p><u>Prepared by</u></p> <p><i>S. Suleet</i> Lab-In Charge <i>Miss S.D. Suleet</i></p>	<p><u>Checked by:</u></p> <p><i>A.G. Jadhav</i> (Prof. A.G. Jadhav) H.O.D</p>	<p><u>Approved by</u></p> <p><i>J.B. Sanap</i> (Dr..J.B.Sanap) Principal</p>
--	--	---