



**GAWANDE COLLEGE OF
PHARMACY.**

Sakharkherda, Lavala Road, Tq.

Sindkhed Raja

Dist.Buldhana



Department:Pharmacognosy

Name of the instrument:Incubator

S.O.P. No.: GCP/PH/030

Stock Register No.:PH0021/08

Effective date:20/02/15

Review date:26/02/15

Procedure:

1. The temperature of incubator is maintained at desired level by thermostat.
2. The incubator is properly ventilated by percolated shape.
3. They are provided with double door.
4. Enter one mode of glass so the contents of incubator may be viewed without administrating outside air.
5. Incubator is used for maintaining constant temperature suitable for growth of various microorganisms

◆ APPLICATION

It is used to carry out growth of microorganism at suitable temperature.

(37°C for bacteria and 27°C for fungi.)



Incubator

Prepared by

S. Seelak

Lab-In Charge

M/SS S.O.D. Solat

Checked by:

A.G. Jadhav

**(Prof.
A.G. Jadhav)**
H.O.D

Approved by

J.B. Sanap

(Dr..J.B.Sanap)
Principal